

Photon Interference

Single Photon and Bi-Photon

K. Oskar Negron

Department of Physics, CSU Fullerton

May 22, 2017

Outline

Photon Interference

K. Oskar
Negron

Introduction

Tools

Single Photon
Interference

Set-Up

Results

Biphoton
Interference

1 Introduction

■ Tools

2 Single Photon Interference

■ Set-Up

■ Results

3 Biphoton Interference

Wave-particle duality

Photon Interference

K. Oskar
Negron

Introduction

Tools

Single Photon Interference

Set-Up
Results

Biphoton Interference

- Light in the form of waves, in or out of phase, either constructively or destructively interfere with one another
- Incident discrete wave packets (photons) with sufficient energy may eject electrons
 - Photoelectric effect
- The fact that light can behave in either of these manners is referred to as Wave Particle Duality.

Outline

Photon Interference

K. Oskar
Negron

Introduction

Tools

Single Photon
Interference

Set-Up
Results

Biphoton
Interference

1 Introduction

■ Tools

2 Single Photon Interference

- Set-Up
- Results

3 Biphoton Interference

Photon Interference

K. Oskar
Negron

Introduction

Tools

Single Photon
Interference

Set-Up
Results

Biphoton
Interference

- The Single Photon Interference Experiment demonstrates that a photon may interfere with itself when it may reach a detector by means of either of 2 indistinguishable paths.
- Feynman provides three Quantum Mechanical Conditions by which interference may occur

3 Quantum Mechanical Conditions Interference

Photon Interference

K. Oskar Negrón

Introduction

Tools

Single Photon Interference

Set-Up

Results

Biphoton Interference

- The probability P of an outcome from a particle interacting with an apparatus is the square of the absolute value of its complex probability amplitude ϕ
 - $P = |\phi|^2$
- If the same outcome can occur in indistinguishable ways, the probability amplitude is the sum of probability amplitudes for each way separately:
 - $P = |\phi_1 + \phi_2|^2$
- If an experiment determines which way the outcome occurred, then the probability is the sum of each alternative
 - $P = P_1 + P_2 = |\phi_1|^2 + |\phi_2|^2$

Can it be explained classically?

Photon Interference

K. Oskar Negrón

Introduction

Tools

Single Photon Interference

Set-Up
Results

Biphoton Interference




Figure: Beam splitting in the classical sense.

- ϵ_2 and ϵ_3 are the amplitude of the reflected and transmitted beams,
- r and t are the complex reflectance and transmittance of the beam-splitter

Can it be explained classically?

Photon Interference

K. Oskar
Negron

Introduction
Tools
Single Photon
Interference
Set-Up
Results
Biphoton
Interference

- $\varepsilon_2 = r\varepsilon_1$ and $\varepsilon_3 = t\varepsilon_1$
- Loss-less beam-splitters
 - $|\varepsilon_1|^2 = |\varepsilon_2|^2 + |\varepsilon_3|^2$
 - $|r|^2 + |t|^2 = 1$
- If treat the beam-splitter quantum mechanically (replace ε_i with annihilation operators \hat{a}_i)

$$\begin{pmatrix} \hat{a}_2 \\ \hat{a}_3 \end{pmatrix} = \begin{pmatrix} r\hat{a}_1 \\ t\hat{a}_1 \end{pmatrix} \quad (1)$$

Can it be explained classically?

Photon Interference

K. Oskar
Negron

Introduction

Tools

Single Photon
Interference

Set-Up
Results

Biphoton
Interference

- Fails the commutation relationships



$$[\hat{a}_2, \hat{a}_2^\dagger] = |r|^2 [\hat{a}_1, \hat{a}_1^\dagger] = |r|^2$$

$$[\hat{a}_3, \hat{a}_3^\dagger] = |t|^2 [\hat{a}_1, \hat{a}_1^\dagger] = |t|^2$$

$$[\hat{a}_2, \hat{a}_3^\dagger] = rt \neq 0$$

What went wrong!?

Correct Interpretation

Photon Interference

K. Oskar
Negron

Introduction

Tools

Single Photon
Interference

Set-Up
Results

Biphoton
Interference




Figure: Beam splitting quantum-mechanically

■

$$\begin{pmatrix} \hat{a}_2 \\ \hat{a}_3 \end{pmatrix} = \begin{pmatrix} r\hat{a}_1 + t'\hat{a}_0 \\ t\hat{a}_1 + r'\hat{a}_0 \end{pmatrix} \quad (2)$$

Mach-Zehnder interferometer

Photon Interference

K. Oskar Negrón

Introduction

Tools

Single Photon Interference

Set-Up

Results

Biphoton Interference




Figure: MZI. BS_1 and BS_2 are beam-splitters, θ is the phase difference, M_1 and M_2 are mirrors, and D_1 and D_2 are detectors.

Mach-Zehnder interferometer

Photon Interference

K. Oskar
Negron

Introduction

Tools

Single Photon
Interference

Set-Up
Results

Biphoton
Interference

- Consider the state $|0\rangle|1\rangle$ entering the MZI and interacting with the first beam splitter (BS_1)



$$|0\rangle|1\rangle \rightarrow t|0\rangle|1\rangle + r|1\rangle|0\rangle.$$

- Once this reaches the second beam-splitter (BS_2) there will/might be a phase length θ_1 and θ_2

$$|0\rangle|1\rangle \rightarrow e^{i\theta_1}(r|0\rangle|1\rangle + t|1\rangle|0\rangle)$$

$$|1\rangle|0\rangle \rightarrow e^{i\theta_2}(t|0\rangle|1\rangle + r|1\rangle|0\rangle)$$

Mach-Zehnder interferometer

Photon Interference

K. Oskar
Negron

Introduction

Tools

Single Photon
Interference

Set-Up
Results

Biphoton
Interference

- the final state is given by,

$$t(e^{i\theta_1}(r|0\rangle|1\rangle + t|1\rangle|0\rangle)) + r(e^{i\theta_2}(t|0\rangle|1\rangle + r|1\rangle|0\rangle))$$

or

$$rt(e^{i\theta_1} + e^{i\theta_2})|0\rangle|1\rangle + (tte^{i\theta_1} + rre^{i\theta_2})|1\rangle|0\rangle.$$

■

$$\begin{aligned} P_{10} &= (e^{i\theta_1} tr + e^{i\theta_2} rt)(\overline{e^{i\theta_1} tr + e^{i\theta_2} rt}) \\ &= 2|t|^2|r|^2 + 2|r|^2|t|^2 \frac{e^{i(\theta_1-\theta_2)} + e^{-i(\theta_1-\theta_2)}}{2} \\ &= 2|t|^2|r|^2 + 2|r|^2|t|^2 \cos(\theta_1 - \theta_2) \end{aligned}$$

Mach-Zehnder interferometer

Photon Interference

K. Oskar
Negron

Introduction

Tools

Single Photon
Interference

Set-Up
Results

Biphoton
Interference



$$P_{10} = 2|t|^2|r|^2 + 2|r|^2|t|^2 \cos \theta$$

- $\theta = \theta_1 - \theta_2$ is the phase difference
- 50:50 beamsplitter have $|r| = |t| = \frac{1}{\sqrt{2}}$



$$P_{10} = \frac{1}{2}(1 + \cos \theta).$$



$$P_{01} = 1 - P_{10} = \frac{1}{2}(1 - \cos \theta)$$



Single Photons?

Photon Interference

K. Oskar
Negron

Introduction


Tools

Single Photon
Interference

Set-Up
Results

Biphoton
Interference

■ Spontaneous parametric down-conversion



Outline

Photon Interference

K. Oskar
Negron

Introduction

Tools

Single Photon
Interference

Set-Up

Results

Biphoton
Interference

1 Introduction

■ Tools

2 Single Photon Interference

■ Set-Up

■ Results

3 Biphoton Interference

Experimental Set-up

Photon Interference

K. Oskar Negron

Introduction

Tools

Single Photon Interference

Set-Up

Results

Biphoton Interference




Figure: Using an arc to determine the locations of the detectors D1 and D2

Experimental Set-up

Photon Interference

K. Oskar
Negron

Introduction

Tools

Single Photon
Interference

Set-Up
Results

Biphoton
Interference

- First, we must assure that we are acquiring coincidence counts in our two detectors that is, they are equidistant from the BBO (Barium Borate) crystal.
- A 405 nm laser is carefully placed to go through the BBO crystal and by means of parametric down conversion, photon pairs are formed 3° apart to be detected in $D1$ and $D2$.
- With the use of a LabView program, our detectors let us know if there has been a photon pair detected in both detectors simultaneously.

Key Procedures

Photon Interference

K. Oskar
Negron

Introduction

Tools

Single Photon
Interference

Set-Up
Results

Biphoton
Interference

- Assure all pieces are the exact same height
- Assure the beam remains that height during its path entirety. (Iris iteration method with EVERY piece)
- Best to have beam remain parallel to holes in breadboard
- Use plumb-bob to accurately place pieces and align beams
- Alignment is done using the red beam to mimic the 405nm beams path
- Optimize coincidence counts

Iris

Photon Interference

K. Oskar
Negron

Introduction


Tools

Single Photon
Interference

Set-Up

Results

Biphoton
Interference



MZI initial set-up

Photon Interference

K. Oskar
Negron

Introduction

Tools

Single Photon Interference

Set-Up Results

Biphoton Interference

- Placing the first beam splitter, assure the incident beam is reflected down its path of incidence and reflected perpendicular relative to incidence
- Assure reflected path is parallel to holes in table using iris iteration
- Place a mirror to reflect the initially transmitted beam parallel to the holes (iterate)
- The first reflected beam hits a mirror on top of a translational stage with a piezo electric to drive the stage very slightly forward, pushing the mirror slightly forward with a small voltage applied.

MZI initial set-up

Photon Interference

K. Oskar Negrón

Introduction

Tools

Single Photon Interference

Set-Up

Results

Biphoton Interference





Figure: Piezo Driver in the Translational Stage

MZI final piece

Photon Interference

K. Oskar
Negron

Introduction

Tools

Single Photon
Interference

Set-Up
Results

Biphoton
Interference

- Place the final beam splitter where the two beams intersect in space.
- Similar to the previous iris iteration process to assure the beams are parallel with the holes and at the correct height
- With the first reflected path blocked, we can align the beam exiting the beam splitter
- Then blocking the first transmitted path, align the translational stage with mirror to superimpose the beam previously aligned and iterate.

MZI final piece

Photon Interference

K. Oskar
Negron

Introduction


Tools

Single Photon
Interference

Set-Up

Results

Biphoton
Interference



MZI final piece

Photon Interference

K. Oskar
Negron

Introduction

Tools

Single Photon
Interference

Set-Up

Results

Biphoton
Interference



MZI with white light

Photon Interference

K. Oskar
Negron

Introduction

Tools

Single Photon
Interference

Set-Up

Results

Biphoton
Interference




Figure: Both the bulb and Fiber Optic Cable are set up

MZI with white light

Photon Interference

K. Oskar
Negron

Introduction

Tools

Single Photon
Interference

Set-Up
Results

Biphoton
Interference

- Determine if the interferometer arms are the exact same length
- Placing a white light source in front of the interferometer and a spectrometer with fiber optic cable to receive the light at the end we can view the white light spectrum.
- If the arm lengths are very near identical, interference will begin.
- Translational stage can be slightly adjusted to an optimal interference of white light or just before that point.

MZI with white light

Photon Interference

K. Oskar Negron

Introduction
Tools
Single Photon Interference
Set-Up
Results
Biphoton Interference




Figure: Some wavelengths interfering

Results

Photon Interference

K. Oskar
Negron

Introduction


Tools

Single Photon
Interference

Set-Up

Results

Biphoton
Interference




Results

Photon Interference

K. Oskar
Negron

Introduction
Tools
Single Photon
Interference
Set-Up
Results
Biphoton
Interference



Biphoton Interference MZI

Photon Interference

K. Oskar Negrón

Introduction


Tools

Single Photon Interference

Set-Up

Results

Biphoton Interference



Biphoton Cases

Photon
Interference

K. Oskar
Negron

Introduction


Tools

Single Photon
Interference


Set-Up
Results

Biphoton
Interference

Case 1: Both of the photons will go through one of the arms MZI.



Case 2: Both of the photons will go through the other arm (different from Case 1) of the MZI.



Biphoton Cases

Photon
Interference

K. Oskar
Negron

Introduction


Tools

Single Photon
Interference


Set-Up
Results

Biphoton
Interference

Case 3: One photon will go through one of the hands of the MZI and the other photon through the other hand of the MZI.



Case 4: Same as Case 3, but the photons pass through the other hand of the MZI.




Complete set-up

Photon Interference

K. Oskar Negrón

Introduction
Tools
Single Photon Interference
Set-Up
Results
Biphoton Interference




Complete set-up

Photon Interference

K. Oskar Negrón

Introduction
Tools
Single Photon Interference
Set-Up
Results
Biphoton Interference



Results

Photon Interference

K. Oskar Negrón

Introduction
Tools
Single Photon Interference
Set-Up
Results
Biphoton Interference

